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
Instructions for use ENGLISH

N1


REF	Désignation / Omschrijving Bezeichnung / Designation	Mesures / Maten / Maßnahmen / Measurements	REF	Désignation / Omschrijving Bezeichnung / Designation
RFB09115	ROLLBELT with 2 side straps S/M	60-110cm	RFP22000	PELVIC
RFB09140	ROLLBELT with 2 side straps L	85-135cm	RFE24000	Lateral extender 55cm
RFB09115A	ROLLBELT with 2 side straps S/M + PELVIC amovible	60-110cm	RFE24500	Lateral long extender 75cm
RFB09140A	ROLLBELT with 2 side straps L + PELVIC amovible	85-135cm	RFB10200	ROLL extender 20cm
RFB09115P	ROLLBELT with 2 side straps S/M + PELVIC	60-110cm	RFB10300	ROLL extender 40cm
RFB09140P	ROLLBELT with 2 side straps L + PELVIC	85-135cm	RFW16100	WRIST HOLDERS 14-22cm
RFB09315	ROLLBELT Comfort S/M with 2 side straps	60-110cm	RFY21000	SHOULDER BELT S/M
RFB09340	ROLLBELT Comfort L with 2 side straps	85-135cm	RFY21500	SHOULDER BELT L
RFB09315P	ROLLBELT Comfort S/M with 2 side straps + PELVIC	60-110cm	RFK01111	MAGNET KEY (M)
RFB09340P	ROLLBELT Comfort L with 2 side straps + PELVIC	85-135cm	RFK01222	MAGNET BUTTON + PIN (M)
RFB10115	ROLLBELT S/M	60-110cm	RFK01201	MAGNET PIN
RFB10140	ROLLBELT L	85-135cm	RFK01202	MAGNET BUTTON
			RFK00100	STEEL KEY (RL)
			RFK00433	STEEL BUTTON + PIN (RL)
			RFK00300	STEEL PIN
			RFK00400	STEEL BUTTON
			REF +M=	BELT + MAGNET LOCKS
			REF +RL=	BELT + STEEL LOCKS

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RENOL products are safety, immobilization, restraint products

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
Indications

1. RENOL products are used for medical purposes to secure an agitated patient in bed, on a stretcher, to avoid an accident, to allow the patient to follow his medical treatment.
2. RENOL products are only used under medical prescription.
3. The products must be used by trained staff following the RENOL's instructions for safe use and according to the establishment's safety rules.
- Contraindications**
4. Ensure that RENOL products do not come into contact with injured skin, do not obstruct breathing nor impede blood pressure.
- Precautions of use**
5. To avoid any accident and in order to ensure the person's comfort and safety, it is essential to read the instruction leaflets carefully and understand them before operating belts. All users follow the instructions and procedures, patients' safety depends on using belts correctly.
6. Complications, accidents related to the restraint in case of misuse, during unsuitable, unsupervised use: bedsores, pneumonia or deep vein thrombosis, psychological stress (agitation, delirium) and negative impact on cognitive abilities, consequence fatal (following severe trauma or suffocation).
7. Ensure that a RENOL product is not damaged before use, all products must be clean and in good condition.
8. We recommend monitoring the immobilised patient frequently; restless and combative patients must be under constant observation.
9. Always use the correct size: too small or too large may compromise safety or comfort. Measure the patient's waist size, wrist and ankle, select the appropriate RENOL product.
10. Remove all sharp objects from the patient (jewellery, watches, glasses, etc.).
11. Only use products RENOL on suitable support. Never use RENOL products as a safety belt in vehicles.
12. When supporting or immobilising the patient in a bed, don't forget the patient's ability to turn over.
13. Connect RENOL products to a solid part of the bed or stretcher. The straps must be secured; they cannot slide or change position; the fastening point of the strap on the bed must not damage the product (sharp edge, etc.).
14. The straps are placed on a part of the bed that can move to match possible movements of the bed base; this also applies to stretchers, wheelchairs and chairs where certain parts can be modified or tipped. Check the straps every time the position changes. The abdominal strap attached to the bed frame should compress the mattress.
15. When using the abdominal belt, always use the pelvic belt RFP22000 and lateral belts RFE24000 unless otherwise specified by the prescribing doctor.
16. Ensure that the RENOL key, part number: RFK00100 (RL) and RFK01111 (M) is always available. Use only RENOL locks.
17. In an emergency, if the key is not available or problem with lock system, cut the strap close to the frame.
18. Before use, check operation of the magnet (M) lock system, check correct operation of the RFK01202 Button with the key RFK01111 and pin RFK01201. We recommend replacing the Button RFK01202 after 2 years. The lifetime may be shorter or longer, depending on use.
19. Before use, check operation of the mechanical lock system (RL), check correct operation of the RFK00400 Button with the key RFK00100 and long pin RFK00300. We recommend replacing the Button RFK00400 every 3 years. The lifetime may be shorter or longer, depending on use.
20. The magnetic key (M) must be placed 10 cm minimum from a pacemaker or any other medical system that could be disturbed by a magnet.
21. The magnetic lock system could be activated with a magnet other than that delivered by RENOL.
- Maintenance and materials**
22. Follow the washing instructions marked on the products.
23. Close Velcro fastenings before washing and drying to prevent them wearing.
24. Detergents: avoid using chlorine.
25. Temperature: maximum temperature 95°C.
26. Drying: delicate cycle at low temperature, about 5% shrinkage possible.
27. Remove all locking systems before washing and drying.
28. Never immerse the magnetic (M) and mechanical closures (RL) in a liquid.
29. Material: 100% Terylene (polyester).
30. The raw materials are treated with fire retardant and anti bacterial.
31. 2 years warranty, this warranty does not cover damage resulting from normal wear and tear.

All our products are CE-approved - Class I medical devices.
Any serious incident occurring in connection with this device must be notified to the manufacturer and to the competent authority of the country where the serious incident took place.

For more information, contact your support specialist, your distributor or www.renolcare.com